

# Standardization of Biomarker Measurements Across ADRCs




Kristen A. Russ, PhD  
NCRAD Biomarker Assay Lab Director  
Assistant Research Professor, MMGE, IUSM

Standardization and quality of biomarker measurements starts with standardizing specimen collection and handling procedures

# NCRAD Protocol Development and Standardization

- Manual of Procedures
  - Includes descriptive and pictorial instructions
  - Training videos available
- Sites have in person or virtual training
  - NCRAD staff walk through the collection and shipping procedures with the team

## Plasma and Buffy Coat Preparation (10ml Lavender-Top Tube x2)



# NCRAD Protocol Development and Standardization

- Sample Verification
  - Provides verification of sample characteristics against the information provided by the site for participants with buffy coats provided
- Blinded Samples
  - All samples provided to labs (including the Biomarker Assay Lab) are received blinded.
- Majority of collection materials are provided in collection kits
- NCRAD receives documentation from sites if samples have non-conformance

# Standardization of biomarker measurements and data handling within the NCRAD Biomarker Assay Laboratory (BAL)

# Specimen Preparation

- BAL technicians are blinded to any information not necessary to track samples as they move through the analysis procedures
  - Site ID, collection date, kit number etc. are not available to the BAL technicians
- NCRAD staff handles creating balanced plate designs to enable BAL to remain blinded
  - Plates/runs are balanced on age, biological sex, and longitudinal samples
  - Samples are provided to BAL in a balanced order

# Specimen Analysis in the NCRAD BAL

## BAL Laboratory Equipment

- Automated liquid handling
  - Tecan Fluent (2)
    - Utilized to prepare Quanterix HD-X plates
      - Lower variability when performing assay preparations compared to onboard processes
    - Regularly monitored using a colorimetric method to ensure consistent and precise readings across all pipettors
    - Regular preventative maintenance as recommended by the manufacturer



# Specimen Analysis in the NCRAD BAL

## BAL Laboratory Equipment

- Immunoassay platforms
  - Quanterix HD-X (2)
    - For use with plasma
  - Fujirebio Lumipulse G1200
    - For use with plasma or CSF
  - Alamar NULISA
    - For use with plasma or CSF





# Assay Selection and Characterization

- Qualification Studies

- Utilizes local sample collection

- 5 samples spanning high, medium, and low concentrations

- 2 days to capture intra and inter-day variability

- Precision
    - Dilutional Linearity
    - Parallelism

- Passing Criteria

- Must have average total CV ~10%
    - Performance must be consistent with the characterization provided by the manufacturer

- Comparability Study

- If assessing new assay for biomarkers currently offered through the BAL, ~100 samples with amyloid or tau PET data are analyzed on the new assay and compared to existing data


# Producing consistent and reliable measurements

Bridging and control measures are key to consistency of data across time

- **Control measures**
  - Identify assay performance outside established quality guidelines
  - Confirms validity of assay results
- **Bridging**
  - Allows comparison of data between studies
  - Necessary to fully harmonize laboratories

# Control Measures

- Assay monitoring within studies and over time at NCRAD
  - Kit QC low and high samples
    - Provided by manufacturer of kits
    - Allows BAL to determine that the assay is working properly from plate to plate
  - Plasma reference pools
    - BAL produces plasma reference pool controls from our local plasma collections
    - High, medium, and low controls
    - Allows monitoring of assay performance within the batch and between batches




Plasma reference pool data points outside  $\pm 2SD$  from the mean are investigated to determine validity of the assay plate/run.

# Bridging

Bridging involves detecting and correcting for lot-to-lot or instrument-to-instrument variability

- Enabled by a local collection of plasma from healthy controls and participants with Alzheimer's disease
- Quanterix HD-X Assays
  - 31 samples run on each lot or instrument
  - CV calculated between the lots
    - Under 10% CV considered acceptable
    - Over 10% requires bridging
  - Also assessed after preventative maintenance or instrument repairs
- Fujirebio Lumipulse Assays
  - 31 samples run between changes for immunoreaction cartridge lots
  - 15 samples run between reagent lot changes
  - CV calculated between the lots
    - Under 10% CV considered acceptable
    - Over 10% requires bridging
  - No bridging required to date
  - Also assessed after preventative maintenance or instrument repairs

# Specimen Analysis in the NCRAD BAL



- Highly SOP driven
- Documentation
  - Assay performance monitoring
  - Sample quality monitoring
  - Process deviations logs
- Standardized data handling
  - Many QC checks
    - QC checks built into data templates
    - QC checks performed at various levels prior to data return
      - Technicians, Coordinator, Lab Director
  - Developing new web app to further standardize and remove human element from data handling

# Specimen Analysis in the NCRAD BAL

**BAL Data Display**

Applications

- Data View
- Data Analysis
- Barcode Grid
- Manifest Status
- Return Documents
- Trackers
- Internal Login

Filters

Assay

Plate

File

Configuration Date

2024-01-10 to 2024-01-10

Save

Users

Config Number

Save Configuration

Select Manifest

☐ Duplicate Barcode Scan?

Scan Type

Initial

Number of Rows

1

Number of Columns

12

☐ Enable/Disable Cell Customization

Save Current State Load Saved State Revert to Default

Barcode Grid Input

A1	A2	A3	A4	A5	A6	A7	A8	A9	A10	A11	A12
B1	B2	B3	B4	B5	B6	B7	B8	B9	B10	B11	B12
C1	C2	C3	C4	C5	C6	C7	C8	C9	C10	C11	C12
D1	D2	D3	D4	D5	D6	D7	D8	D9	D10	D11	D12
E1	E2	E3	E4	E5	E6	E7	E8	E9	E10	E11	E12
F1	F2	F3	F4	F5	F6	F7	F8	F9	F10	F11	F12
G1	G2	G3	G4	G5	G6	G7	G8	G9	G10	G11	G12
H1	H2	H3	H4	H5	H6	H7	H8	H9	H10	H11	H12

Download ID Barcodes

- New web app will further standardize data handing and QC
  - Limits data manipulation
  - Simplifies and decreases time for data analysis and QC
  - Streamlines monitoring of assay performance

